Ionomolar

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Ionomolar is a radiopaque light curing **compomer filling material**. It contains a radiopaque glass of fluorosilicates in a matrix of polymerisable monomers and other light curing polymers.

Ionomolar combines the advantages of glassionomer cements - good adhesion to tooth material, releasing of fluoride and biocompatibility - with those of light curing composites, namely mechanical stability, esthetic and color stability. After light curing the filling material has a high gloss polishability.

Indications

- restoration of teeth of the first dentition
- restoration of class III cavities
- restoration of class V cavities
- restoration of cervical abrasions and erosions

Side Effects

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection the cavity floor in cases of deep excavations should be covered with a thin layer of calciumhydroxide material.

Incompatibility with Other Materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the compomer. Neither store the composite material in proximity of eugenol containing products, nor let the compomer allow coming into contact with materials containing eugenol.

Color matching

Clean the tooth with a fluoride-free polishing paste (e.g. flour of pumice with water) prior to preparation and color matching. Ascertain the tooth shade while teeth are still moist and select the appropriate Ionomolar shade. Ionomolar is shaded according to VITA®-Shades.

Preparation

In all cavity classes cavity is prepared to the minimum required for caries removal.

Cavity cleanliness is essential for the adhesiveness and durability of the compomer filling. In cases where no cavity preparation is necessary clean the tooth surface with flour of pumice and rotating rubber instruments.

Wash thoroughly with air/water spray and dry with oil free air. **Do not desiccate** the dentine.

For improving the results enamel and dentine etching (**15 sec**) is recommended (e.g. with Dentamed Etching Gel). After etching rinse with water for **15 sec** and dry for **10 sec** in an oilfree stream of air. **Do not desiccate** the dentine.

Avoid contamination with saliva once cleaning has been carried out. Therefore isolating of the working field with cotton rolls or a rubber dam, depending on the clinical situation.

Pulp protection

Cavity floor of deep excavations should be covered with a thin layer of calciumhydroxide material.

Application

After cavity preparation as above described and application of a suitable light cure dentin/enamel adhesive and bonding system (e.g. Dentamed P&B II) Ionomolar can be applied in the selected shade.

Ionomolar is easy to model by using suitable instruments. Transparent matrix strips may be used.

Curing of Ionomolar

For an optimum result apply **Ionomolar** in **layers of max. 2 mm**. Light cure each increment separately.

By using a polymerization unit with a **light intensity of at least 1000mW/cm²** and a **wavelength range of 400 to 500 nm** (blue light) **cure each increment as follows:**

Lighter shades (e.g. A1, A2, B2) **20 seconds**

Darker shades (e.g. A3.5, B3) **30 seconds**

Finishing

Contour with fine diamonds, stones or burs. Polish to high gloss with discs or rubber points. Interproximal finishing is accomplished by fine grit finishing strips.

Storage

Store at temperatures below 25 °C. Keep in a dry place. Store in the original seal-pack until use. Once the seal-pack has been opened, Ionomolar may be stored for a maximum of 2 months.

Additional Notes

- Do not use any resin or other compomer materials to adjust viscosity of Ionomolar material.
- The product is **not suitable** for direct or indirect overcapping of the pulpa and for occlusal stress-bearing restorations of the second dentition.
- Ionomolar contains methacrylates, that may be irritating to the eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical device.
- For susceptible patients sensitization to the product by contact with skin or mucous membranes can not be excluded. After accidental contact, wash immediately with soap and water or rinse off with plenty of water.
- If patient reacts allergic, don't use Ionomolar any more. If fillings of Ionomolar has to be removed and it is surely known, that the patient reacts allergic to the product, use of protecting rubber dam is absolutely advised to prevent occlusions of particles in the mucous membrane. Prophylactic removal without protection by help of rubber dam is therefore strictly contraindicated.
- APPLIC.-Tips are for single use only.
- VITA® is a registered trade mark of the VITA-Zahnfabrik, Bad Säckingen, Germany.

Warranty

P.L. Superior Dental Materials GmbH warrants this product will be free from defects in material and manufacture. P.L. Superior Dental Materials GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and P.L. Superior Dental Materials GmbH's sole obligation shall be repair or replacement of the P.L. Superior product.

Limitation of Liability

Except where prohibited by law, P.L. Superior Dental Materials GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep away from children!**For dental use only!**

Distributor:

Dr. Ihde Dental GmbH
Erfurter Strasse 19
D-85386 Eching
Tel. +49 89 319 761-0
info@ihde-dental.de