Deutsch:

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English:

Instructions for use – Please read carefully

IHDE® GIZ® fil capsules

IHDE® GIZ® fil capsules is a glass polyalkenoate filling material in capsule form, which is used, for example, for definitive restorations of cavity classes III and V, milk tooth cavities, tunnel preparations, root caries and liners. Over a longer period of time, small quantities of fluoride ions can get into the solution as water enters. They change the apatite of the surrounding tooth structure into fluoride apatite and make it more resistant to caries. IHDE® GIZ® fil capsules bonds well to dentin and enamel because of the reaction of the carboxylate groups of the polyacrylic acid and the calcium of the hydroxylapatite. Considerably more healthy tooth material can be retained as a result of this direct chemical bond because only flat undercuts are necessary.

Indications

- Fillings for cavities of Black's class I in milk teeth
- Filling of cavities of Black's class III and V including root erosions and wedge-shaped defects
- Stump build-ups
- Fissure sealing
- Treatment of caries at approximal contact points (tunnel preparation)
- Semi-permanent restorations of cavity classes I and II

Recommended use

1. Preparation of the cavity

As much healthy tooth material as possible should be preserved because only flat undercuts are required for a better bond. The edge should not be thinner than 0.5mm.

2. Drying the area for treatment

The area to be treated should be covered with a matrix band, if possible.

3. Pulp protection

Rinse the prepared cavity thoroughly with water and dry in an oil-free current. Calcium hydroxide preparations should be used to protect the pulp in very deep cavities.

4. Pretreatment of the dentin

Apply a tooth cleaner to the recently prepared dentin surfaces to remove the lubrication layer. Follow the manufacturer's instructions. Avoid dehydration when drying with oil-free air!

5. Activation of von IHDE® GIZ® fil capsules

Activate the capsule in the application forceps by squeezing once and agitate for 10 seconds in a suitable capsule mixing appliance at around 4000 oscillations per minute.

6. Application of IHDE® GIZ® fil capsules

Apply the mixed cement to the cavity directly from the capsule by squeezing the application forceps several times. Pack the cement down firmly into the cavity. If necessary, shape the surface with an appropriate matrix. Use a protective varnish for the first 24 hours after application to stop water from getting into the filling. IHDE® GIZ® fil capsules can be worked for 2 minutes (including mixing). Do not work the material once the hardening process has begun. Hardening takes 3 minutes.

7. Finishing

After 3 minutes the matrix can be removed and this area, which is now unprotected, can be coated with a tooth varnish. Finishing should begin after a further 6 minutes. Excesses should be removed with Vaseline-coated corundum grinders or sandpaper disc and the filling polished. However, polishing is only indicated on the edge of the filling.

Please note

- Do not use water spray!
- Protect areas of dentin around the pulp
- Dry the area to be treated.
- Caution: Do not allow uncured cement to come into contact with the eyes. In Case of emergency, rinse immediately with plenty of water.

Contra-indications

Use of IHDE® GIZ® fil capsules is contra indicated:

- In cases of direct application to dentin around the pulp or to exposed pulp
- if a dry working area or the prescribed application technique is not possible
- In fillings with high chewing load
- In cases of proven allergies to the ingredients of IHDE® GIZ® fil capsules

Side effects

Serious irritation of the pulp may occur with use of glass polyalkenoate cements. Side effects of this kind have not been experienced in IHDE® GIZ® fil capsules.

Composition

IHDE® GIZ® fil capsules powder contains an ionomer glass and dried polyacrylic acid. IHDE® GIZ® fil capsules liquid contains aqueous polyacrylic acid solution, polycarbonate acid and tartaric acid.

Classification

IHDE® GIZ® fil capsules meets requirements in accordance with ISO 9917-1:2003.

Storage

Store in a dry place and protect against moisture. Do not store over 25°C! Do not use after the expiration date.

An open packet may only be stored for a maximum of 6 months (depending on the expiration date of the unopened packet), if stored as prescribed.

Keep away from children. For dental use only.

Packaging units

50 capsules @ 306mg powder and 85mg liquid
 Universal A2 
 Yellow B3 

Legends

-  Read instructions
-  Expiration date
-  LOT Charge number
-  Keep in a dry place
-  Store tightly closed keep
-  Temperature range from 15 °C to 25 °C
-  Manufacturer

CE 0481

9999-2001-GA-01

Revision issued 2018-07-05

 **Dr. Ihde Dental GmbH**
 Erfurter Str. 19
 D-85386 Eching/München
 Tel. +49 (0)89 319 761-0
 Fax +49 (0)89 319 761-33
 info@ihde-dental.de
 www.ihde-dental.de