English:

Instructions – Please read carefully

## IHDENT® CemBase® blue (Powder)

CemBase® is a water-mixable glass-ionomer cement used for cavity liners, core build-ups and cementation in prosthetic, conservative and orthodontic dentistry. It is quick setting and fluoride releasing.

**Essential components**

Aluminium silicate fluoride, polyacrylic acid, iron oxide pigments (known)

**Areas of application**
Cementation of prosthetic restorations and orthodontic appliances (e.g. crowns, bridges, inlays/ onlays, bands)
Cavity liners

**Interactions**

None known

**Side-effects**

Severe irritation of the pulp may occur.

**Note:** It cannot be ruled out that CemBase® may cause sensitisation in sensitive patients. CemBase® should no longer be used if these symptoms occur. CemBase® should not be used if there is a suspected allergy. If it is necessary to remove the CemBase® in such cases, a rubber dam should be placed before removal to prevent components of the material infiltrating the oral mucosa. It is essential to apply calcium hydroxide (**REF** 8768) with deep cavities and minimal residual dentine layer (dentine third near the pulp) to prevent possible pulp reaction.

Over dispensing the water may lead to increased solubility of the cement and consequently pulp reaction. Hydrostatic pressure in close fitting crowns can also cause pulp reaction. Excessive over filling of crowns should therefore be avoided.

**Contraindications**

Do not apply directly to exposed pulp. Use is only indicated if it is possible to ensure relative isolation.

**Caution**

Avoid contact with the eyes.

**Application - Preparation**

Carefully clean and dry the enamel, dentine and metal surfaces to ensure an optimal bond of Cem-Base®. Avoid overdrying, e.g. with an air syringe. It is sufficient to dry the surface using a pledget. A moisture film, as is usual with vital teeth, does not impair setting of the cement. Condition the dentine surface for 10 seconds using IHDENT® GIZ® Conditioner (**REF** 2330) to improve the bond strength. Then rinse off the conditioner and dry the cavity. Avoid recontamination with water or saliva after this cleaning stage.

**Pulp protection**

It is essential to apply calcium hydroxide (**REF** 8768) with deep cavities and minimal residual dentine layer (dentine third near the pulp). The use of a conditioner is contraindicated in this area.

**Dispensing**

Minimum 2 scoops of powder and 1 drop of water. The powder and water should be mixed on a mixing pad or clean glass plate. To ensure accurate dispensing, fill the spoon liberally with powder and wipe off the excess on the neck of the bottle when removing the spoon from the bottle. Normal tap water can also be used. Hold the water bottle vertically when dispensing. Adhere strictly to the powder/water ratio for cementation. If a lower viscosity is required, e.g. for cavity liners, this can be attained by increasing the proportion of water accordingly. Mix on a cooled glass plate to compensate for the resulting shorter working time.

**Mixing**

Abrasion-resistant metal cement spatulas can be used for mixing. It is not necessary to use plastic spatulas. The powder and water should be mixed on a mixing pad or glass plate. The powder should be worked into the liquid to produce a viscous consistency. (The required consistency for cementation is attained when the cement only just drops from the spatula.) We recommend mixing portions of powder into the liquid at a time (e.g. in 3 portions) to attain better wetting.

**Working time**

The following times apply for the individual working stages at a room temperature of 23 °C and 50 % humidity:

Mixing time	<b>15 seconds</b>
Working time	<b>2:25 minutes</b>
Setting time	<b>2:25 minutes</b>
Total setting time	<b>4:45 minutes</b>

Higher temperatures shorten the working time and lower temperatures extend the working time. Exceeding the working time will result in a loss of adhesion to the enamel and dentine.

**Application**

Apply a thin layer uniformly to the prepared tooth and crown margins and fit the restoration to its final position without applying pressure. Avoid excessive overfilling. Excess cement should only be removed once the cement has fully set. Protect the filling against water penetration with IHDENT® GIZ® Varnish (**REF** 2440). These cement segments should be only be removed once the cement has fully set.

**Note:** CemBase® adheres to metal instruments and should therefore be washed off using cold water before setting.

**Storage**

The maximum storage temperature is 25 °C. Store dry, tightly close the container after use. Refer to the barcode for batch number and expiry date.

**Keep out of reach of children. Intended only for use by dentists.**

**Pack sizes**  
Pack with 50 g powder, accessories **REF** 8521

<b>Legend</b>	
	Read instructions
	Expiry date
	LOT Charge number
	Keep in a dry place
	Store tightly closed keep
	Temperatures up to 25 °C
	Manufacturer

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IHDENT®, CemBase® and GIZ® are registered trademarks.