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Instructions for Use**Safeprint MV**

Super hydrophilic impression material for the (one step) two phase impression technique as well as one phase impression material

Safeprint MV is very versatile. It can be used:

- as a hydrophilic syringable material for the (one step) two phase impression technique in conjunction with Safeprint P
- as a Correction material in the (two step) putty/wash technique
- as a syringable One Phase Impression Material for crowns, bridges, inlays, partial and complete dentures and denture repairs. It also serves as an impression material for relining, functional impressions, impressions for prosthetic restorations and for impressions using the copper band technique.

Safeprint MV is a light to medium viscosity hydrophilic precision impression material based on addition curing vinyl silicones delivered in Automix-cartridges. It can be extruded directly onto the teeth or the tray. Safeprint MV shows strong hydrophilic properties. Thus improves the wettability and yields an exact reproduction of subgingival areas and other moist surfaces.

Safeprint MV is thixotropic. That means it stays in the position in the tray, but under pressure it flows. It meets the requirements of ISO 4823 Type 3 (light-bodied consistency).

Gun load and extruding

Load the gun with the Automix cartridge and mixing cannula.

For direct extruding onto teeth use "intra oral tips".

Impression taking

Load the auto mixed material onto the tray or directly onto the teeth. Insert the tray not later than **2 minutes** (working time at 23°C or 74°F) after extruding Safeprint MV into the mouth. Maintain position for at least **3 minutes**.

Notes:

Store used cartridge with fixed used mixing cannula.

Do not re-use the cartridge cap. Closing the cartridge by turning on the original cap between use could cause cross contamination of components and thus premature setting.

Casting

The impression should be stored at normal room temperature (max. 25 °C or 77 °F). The impression can be poured **30 minutes** after removal from the mouth. No other time limitation does exist.

Electroplating

All impressions can be copper- or silverplated.

Desinfection

Impressions can be disinfected with standard desinfection solutions.

Storage

Do not store above 25 °C (77 °F). Do not use after expiry date.

Additional notes:

Do not use with standard (condensation reaction type) silicone impression materials.

Do not use Latex gloves, they may interfere with the setting reaction of the impression material.

Warranty

P.L. Superior Dental Materials GmbH warrants this product will be free from defects in material and manufacture. P.L. Superior Dental Materials GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and P.L. Superior Dental Materials GmbH's sole obligation shall be repair or replacement of the P.L. Superior product.

Limitation of Liability

Except where prohibited by law, P.L. Superior Dental Materials GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep away from children!

For dental use only!

Distributor:

Dr. Ihde Dental GmbH
Erfurter Strasse 19
D-85386 Eching
Tel. +49 89 319 761-0
info@ihde-dental.de

Manufacturer:

P.L. Superior Dental Materials GmbH, 22607 Hamburg, Germany