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**Instructions for Use - Please read carefully**

**Drills**

**Instructions for Drill**

For proper use of our products it is recommended that practitioner is being trained. Dr. Ihde Dental AG organizes periodic training courses, which reports on the particularities of each of the components of surgical instruments and their proper use. Instructions for use are available in pdf format at the following web address: [www.simpladent-implant.com](http://www.simpladent-implant.com) No. 31-3  
To open pdf files you need the free Adobe Acrobat Reader.

**Intended use**

The "Surgical Drills" for the implant system are intended to be used for the preparation of the implant bed for inserting dental implant systems.

**Surgical instruments description**

For proper selection of our products are advised to consult the System Application Propectus (SAP). All surgical instruments placed on the market by Dr. Ihde Dental AG are made of stainless steel. Surgical instruments for preparation of the implant site, drills are provided with depth markings which enable the user to know the depth at which the bone is being drilled. The depth indication that drills does not take into account the tip, this length is being proportional to the diameter of the cutter.

**PRESENTATION OF DR. IHDE DENTAL AG SURGICAL INSTRUMENTS**

**Packaging**

Dr. Ihde Dental AG ensures that all surgical instruments follow a manufacturing process, control and extreme cleaning before being packaged. Surgical instruments are packed in blister, bag or box. The product contains identifying label with the batch number, size and model of product.

**Instrument supply state**

Surgical instruments are not sterile. Before using, product must be sterilized, at customer sites. For the correct sterilization of the product you should follow the following steps:

- Removing the product from the packaging supplied. This package does not allow the proper sterilization of the product contained.
  - Place the product in e.g. the Surgical kit or any other suitable packaging or container. The drills are to be used for the recommended sterilization parameters.
- We recommend steam sterilize (autoclave) using a pre-vacuum cycle. Please observe the characteristics of equipment you have. The parameters recommended by ANSI / AAMI ST79 are:
- Temperature: 270°F / 132°C
  - Exposure time: 4 minutes.
  - Drying time: 20 minutes.

**Surgical Kits**

Dr. Ihde Dental AG supplies different models of kits that allow a rational organization, and proper storage of surgical instruments. Surgical kits or trays are manufactured from high performance thermoplastic materials which allow undergo the process of steam sterilization in an autoclave. For proper sterilization is advised to follow the parameters listed in the previous section.

**Contradictions**

None Known.

**Warnings and precautions**

Proper irrigation is required while using surgical instruments in order to prevent overheating of the bone which could cause failure of the implant.

Surgical instruments are part of an overall design concept, having adequate use with products manufactured by Dr. Ihde Dental AG, therefore, the combination of Dr. Ihde Dental AG's products with others outside the mark can cause no optimal operation of these products.

Surgical instruments equipped contra-angle attachment, as ISO 1797, should not be employed at a torque higher than 70 Ncm. The surgical instruments must be cleaned and sterilized before use on another patient. Products should be used only for the purpose for which they were specified.

Do not use the drills if there is a known or suspected allergy against the drill material (stainless steel). Before proceeding with surgery, it is the user's responsibility to check the packaging, the container and to make sure it is the right product for the patient.


For proper cleaning of surgical instruments the following instructions should be followed:

- Do not use cleaning or disinfection with chlorine or acids. Do not use surface disinfectants on our drills.
- Avoid all products containing an aldehyde, which the ability to attach proteins.
- Remove carefully all post-operative debris left on the surface or inside the surgical instruments. For this work use a nylon bristle brush.
- For the instruments with cutting edges, to avoid as much as possible the contacts with other instruments that may damage the surface and / or the cutting edges.
- Immediately after proceeding to the cleaning process, rinse with distilled water and dry thoroughly.
- Before using a new disinfectant agent, try on one sample drill and make sure that the color-coding does not get dissolved.
- Improper drying can cause occurrence of oxidation points in the products.

As the surgical drills are the subject to wear and use, the responsibility for the periodic renewal of these products remains on the user.

Dr. Ihde Dental AG advises to check the status of the drills before surgery and in dependence of the bone type recommended the replacement after 10 to 20 uses, depending on the quality of the bone. In general, during the oral use of our products the user must take precautions to prevent aspiration by the patient as well as overheating of the bone.

**Legend**

	Read instructions
	Expiration date
	Only use once
	non sterile
	This product is only for sale to dentists, oral surgeons and maxillo-facial surgeons or on their behalf.
	Charge number
	Store in dry place
	Manufacturer
	Catalogue number

**Liability, safety and warranty**

If the protective packaging of product is damaged during delivery, Dr. Ihde Dental AG will replace it free of charge. This will not be the case if there is evidence of product manipulation.

The guarantee will apply to Dr. Ihde Dental AG instruments provided they have been used following the instructions indicated. For a refund of the products that are mentioned in this instruction the guidelines established in our conditions of sale and delivery have to be followed.

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**Dr. Ihde Dental GmbH**
**onewaybiomed GmbH**