Deutsch:

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English:

Instructions for use – Please read carefully

IHDENT® P&B® II

IHDENT® P&B® II is a light curing single component primer & bonding with water base.

IHDENT® P&B® II is easy to use and achieves high adhesive values on dentin, enamel and metal (especially amalgam). The adhesive quality nearly corresponds to that of glass ionomer cement.

IHDENT® P&B® II is different than other bonding systems: Highly biocompatible adapted methacrylic ester creates very good adhesion, significantly higher than common bondings based on the phosphate groups.

IHDENT® P&B® II can however be combined with all common light curing composites. The use with Triomolar® REF 4003 or Enamelite® REF 4050 is recommended.

Constituents

Modified polyacrylic acid, hydrophilic methacrylic monomers, camphor-chinon, stabilizers etc.

Use

Drying

The use of a cofferdam is recommended in the drying. Cavity preparation. According to the rules the preparation for composite cavities can be minimally prepared. The margins should be slightly chamfered (0.5 - 1 mm). Deep cavities should be covered with a Ca(OH) preparation in the pulpal direction.

Etching

Phosphoric acid etching gel should be applied in the area of the enamel and the dentin and etched for 20 seconds. Deciduous teeth and teeth with a high fluoride content are etched for 90-120 seconds. Afterwards rinse thoroughly and dry with oil free air. The enamel now appears chalky.

Note

- The etched surface may not be contaminated prior to the application of IHDENT® P&B® II. Should this however occur, the etching process must be repeated.
- Avoid the contamination of skin and mucous membranes with the etching gel. In case, rinse thoroughly with water.
- IHDENT® P&B® II is worked into the tooth surface during the 60 seconds with a brush. If necessary, apply a second layer in order to assure a good coating of the surface. Then the tooth surface is dried with the air blower until it really appears to be dry. Afterwards it is illuminated with the polymerization lamp for 20 seconds without previously applying composite. Now the composite can be applied as necessary.
- IHDENT® P&B® II cures under lack of oxygen. This occurs in the dentin tubuli or in the pores on the enamel surfaces. IHDENT® P&B® II remains viscous on the locations receiving oxygen.

Warning

IHDENT® P&B® II contains polymerizable monomers which can induce allergic reactions in sensitive persons. In this case, the use is to be discontinued and the filling therapy with composites to be contemplated in this patient. The removal of composite fillings may only be performed with a cofferdam in such cases in order to prevent an inoculation of filling and bonding components in the mucous membrane.

Precautionary statements

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Side effects

If sensibilizations occur, all accessible remains of IHDENT® P&B® II should be carefully removed from the skin and mucous membrane and the material should not be used in this patient.

Shelf life


3 years after production and upon storage at 15-25 °C. Store cool, dry and protected from light.

Keep out of reach of children. Only to be used by the dentist.

Package sizes

Bottle with 10 ml IHDENT® P&B® II REF 5125
Bottle with 3 ml IHDENT® P&B® II REF 5126

Legends

	Read instructions
	Expiration date
	LOT Charge number
	Keep in a dry place
	Store tightly closed keep
	Temperature range from 15 °C to 25 °C
	Manufacturer

CE 0029

Warning

- H315**
Causes skin irritation.
- H319**
Causes serious eye irritation.
- H335**
May cause respiratory irritation.

9999-5125-03

Revision issued 2015-07

Dr. Ihde Dental AG

Dorfplatz 11
CH-8737 Gommiswald/SG
Tel. +41 (0)55 293 23 23
Fax +41 (0)55 293 23 00
contact@implant.com
www.implant.com

Distribution by

Dr. Ihde Dental GmbH

Erfurter Str. 19
D-85386 Eching/München
Tel. +49 (0)89 319 761-0
Fax +49 (0)89 319 761-33
info@ihde-dental.de
www.ihde-dental.de

IHDENT®, P&B®, Triomolar® and Enamelite® are registered trademarks.

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**Dr. Ihde Dental AG**Dorfplatz 11