Deutsch:

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English:

Instructions for use – Please read carefully

IHDENT® P&B® III

IHDENT® P&B® III is a filled, light curing single-component bonding agent for enamel and dentin bonding in conjunction with the acid etch technique.

The most important constituent

IHDENT® P&B® III contains HEMA, dimethacrylates, phosphonic acid acrylate, highly dispersed silicon dioxide, initiators and stabilizers in an alcohol solution.

Indication

- Adhesive for direct, light curing composite and compomer restorations
- Adhesive for light initiated resin cements used in conjunction with indirect all ceramic and composite restorations (venners, inlays, onlays, and full coverage restorations)

Contraindication

Do not use IHDENT® P&B® III if a patient is known to be allergic to any of the materials ingredients or if the stipulated working technique cannot be used. IHDENT® P&B® III is a light-curing adhesive only and should not be used in conjunction with self-curing resin materials. Dual-curing materials must be light initiated.

Side effects

IHDENT® P&B® III may cause a sensitizing reaction in patients with a hypersensitivity to any of its ingredients. In these cases, the material should no longer be used.

Interaction

Materials containing eugenol or clove oil may inhibit the polymerization of IHDENT® P&B® III.

Application

Direct restorations (composites, compomers)

1. Apply phosphoric acid gel to the prepared enamel and then flow the etchant onto the prepared dentin. Use a disposable brush or syringe tip to push the etchant into any restricted areas of the preparation, leave for 15 - 30 seconds on enamel and for 10 - 15 seconds on dentin. Remove all etchant gel with a vigorous water spray for at least 5 seconds. Excess water is to be removed until the surface of dentin should have a uniform, glossy appearance (wet bonding). **Do not overdry dentin!**

2. Saturate enamel and dentin with a generous amount of IHDENT® P&B® III onto all prepared dentin surfaces and agitate gently for at least 10 seconds. Take care that the cavity walls are completely covered. Remove excess material. Cure IHDENT® P&B® III with halogen, plasma arc, or laser polymerization lights. When using high-performance halogen curing light, the polymerization time is 10 seconds. For other halogen curing lights, the polymerization time is 20 seconds. The curing time with plasma arc or laser polymerization lights depends on the intensity and the emitted wavelength range of the corresponding curing light (contact the manufacturer of the curing light, if necessary). **Polymerization times below 5 seconds are not recommended.**

Indirect restorations (ceramics, composites in a suitable cementation system)

1. Prepare the surface of the ceramic or composite restoration according to the instructions of the manufacturer. Apply IHDENT® P&B® III onto the internal aspects of the restoration and disperse the material to a thin layer clean, dry air. Do not light-cure the adhesive at this step and protect the restoration from light.

2. Apply phosphoric acid gel (37 %) to the prepared enamel and then flow the etchant onto the prepared dentin. Use a disposable brush or syringe tip to push the etchant into any restricted areas of the preparation, leave for 15 - 30 seconds on enamel and for 10 - 15 seconds on dentin. Remove all etchant gel with a vigorous water spray for at least 5 seconds. Excess water is to be removed until the surface of dentin should have a uniform, glossy appearance (wet bonding). Do not overdry dentin; moist preparation surfaces should be apparent.

3. Saturate enamel and dentin with a generous amount of IHDENT® P&B® III and agitate gently the adhesive onto all prepared dentin surfaces for at least 10 seconds. Take care that the cavity walls are completely covered. Remove the excessive material. Then cure IHDENT® P&B® III with halogen, plasma arc or laser polymerization lights. When using high-performance halogen curing light, the polymerization time is 10 seconds. For other halogen curing lights, the polymerization time is 20 seconds. The curing time with plasma arc or laser polymerization lights depends on the intensity and the emitted wavelength range of the corresponding light (contact the manufacturer of the curing light, if necessary). **Polymerization times below 5 seconds are not recommended.** Subsequently, the restoration is seated using a light / dual-curing luting cement. The materials are then light-cured from all surfaces for at least 40 seconds each according to manufacturers instructions. Proper output of operator curing light must be maintained at all times.

Precautionary statements

- P260** Do not breathe dust / fume / gas / mist / vapours / spray.
P280 Wear protective gloves / protective clothing / eye protection / face protection.
P301+P330+P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water / shower.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P333+P313 If skin irritation or rash occurs: Get medical advice / attention.

Keep away from children. For use in dentistry only.

Storage

Batch-no. and expiry date see label / packaging. Do not use after the indicated date of expiration. Storage temperature: 2 - 28 °C / 36 - 82 °F. Close bottle immediately after use.

Presentation

3 ml bottle IHDENT® P&B® III

REF 5123

Legends

-  Read instructions
-  Expiration date
-  LOT Charge number
-  Keep in a dry place
-  Store tightly closed keep
-  Temperature range from 2 °C to 28 °C
-  Manufacturer

CE 0029

Danger

H314 Causes severe skin burns and eye damage.

H317 May cause an allergic skin reaction.

9999-5123-03

Revision issued

2015-07

Dr. Ihde Dental AG
 Dorfplatz 11
 CH-8737 Gommiswald/SG
 Tel. +41 (0)55 293 23 23
 Fax +41 (0)55 293 23 00
 contact@implant.com
 www.implant.com

Distribution by

Dr. Ihde Dental GmbH
 Erfurter Str. 19
 D-85386 Eching/München
 Tel. +49 (0)89 319 761-0
 Fax +49 (0)89 319 761-33
 info@ihde-dental.de
 www.ihde-dental.de

IHDENT® and P&B® are registered trademarks.

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