Deutsch:

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English:

Instructions for use – Please read carefully

## IHDENT® GIZ® fil capsules

IHDENT® GIZ® fil is a chemically curing, water-mixable glass-ionomer cement used as a filling material in restorative dentistry, as an orthodontic cement and as a cement for build-ups in prosthodontics. IHDENT® GIZ® fil is a fluoride donor.

### Essential components

Ionomer glass, polyacrylic acid, water additives

**1 capsule contains** 306 mg of powder  
85 mg of liquid

Net volume per capsule: 255 mg and 0.1ml paste respectively.

### Application areas

- Class V restoration in cervical lesions of mechanical or abrasive origin
- Class III restorations, especially in large lesions and in lower incisors
- Restorations of all classes in deciduous teeth
- Carious root lesions
- As a fissure sealant for premolars and molars

### Special advantages

- No undercuts required thanks to chemical adhesion to enamel and dentin
- No discolouration of the restoration
- Reliable marginal seal
- Rapid hardening, allowing rapid loading
- Excellent translucency, good colour adaptation
- Good pulp tolerance, since a near-neutral pH is attained shortly after mixing and no acid overdose is possible

### Adverse reactions

Glass-ionomer cements can be strongly irritating to the dental pulp. No such adverse effects have been reported for IHDENT® GIZ® fil to date. Avoid contact of the paste with skin. If accidental contact with skin occurs, wash skin immediately with soap and plenty of water and consult a physician.

### Contraindications

Use of IHDENT® GIZ® fil is contraindicated for direct application to dentin close to the dental pulp or to the exposed dental pulp itself, for restorations subject to high masticatory loads and in patients known to be allergic to one of the ingredients of IHDENT® GIZ® fil. The placement is contraindicated when a dry working field or the recommended application technique is not possible.

### Application

#### Preparing the cavities

Retain as much healthy tooth substance as possible, as no undercuts are required to improve retention. Margins should not be thinner than 0.5 mm.

#### Pulp protection

Rinse the prepared cavities thoroughly with water and dry with oil-free pressurized air. For very deep cavities, use a calcium hydroxide preparation (REF 8768) as a liner.

#### Preparing the dentin

To improve chemical adhesion, the surface of the tooth should be treated with a conditioner. Use a tip to apply the conditioner (REF 2330) or a different polyacrylic acid solution. Leave in place for 10 seconds, rinse thoroughly with water and dry with pressurized air. Avoid excessive drying.

#### Activating and mixing the capsule

Insert the capsule by letting its groove engage the holder of the applicator pliers. Activate the capsule by firmly closing the pliers once. This will loosen the clip, allowing the powder and liquid to mix. Insert the capsule into the mixing fork of a capsule mixer and mix at approximately 4000-4500 cpm. The mixing time is 10 seconds.

#### Applying IHDENT® GIZ® fil

Insert the capsule into the application pliers. Extrude the mixed material by multiple pressing. Introduce the mixed material directly into the cavity and condense. Shape the surface with an appropriate matrix if required and coat with an ionomer varnish (such as Varnish, REF 2440). The working time for IHDENT® GIZ® fil is 2 minutes including mixing. Do not manipulate the material once the setting process has begun. The setting time is 3 to 3.5 minutes. To remove the empty capsule, press the reset button on the application pliers.

**Caution:** Keep glass ionomer powder, dental cleaner, or ionomer varnish away from the eyes in order to avoid chemical burns caused by the polyacrylic acid. If contact with the eye has occurred, immediately and carefully rinse with copious amounts of water and consult an ophthalmologist.

**Note:** Immediately clean all instruments under water after use.

#### Finishing and polishing

Once the material has hardened, you can remove the matrix. Coat the unprotected aspects with ionomer varnish and start finishing after an additional 6 minutes. Remove excess material using Vaseline-coated abrasive stones or flexible cutting wheels. Perform final polishing with Vaseline-coated polishing stones or wheels.

**Caution:** Do not use any water spray.

Coat the freshly prepared restoration with another layer of ionomer varnish to protect the glass-ionomer cement from water during the first 24 hours. Do not wet-polish until after 24 hours. Once the polishing has been completed, the restoration can once again be protected with varnish to protect the continuing hardening process.

#### Orthodontic use

Once the orthodontic band has been adapted, dry the appropriate tooth and carefully coat the band with IHDENT® GIZ® fil on the inside. Then immediately place the band. You must remove any cement excess during the short viscous-plastic phase of the setting process. The position of the band must not be changed. To remove the band, loosen it and take the band off. Remove residual cement using periodontal instruments of ultrasonic devices.

### Precautionary statements

**P280** Wear protective gloves / protective clothing / eye protection / face protection.

**Keep out of the reach of children. Intended for use by dentists only.**

**In order to remove the capsule from the application instrument, press the knob.**

**In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.**



### Storage / Shelf life

Store tightly closed in a dry place. Maximum storage temperature: 25 °C. Check the label or bar code on the packaging for batch number and expiry date. Do not use after expiry date. Shelf life is 6 month when the package has been opened; provided it has been stored dry and at room temperature.

### Package sizes

Package of 50 capsules  
universal A2 REF 2001  
yellow B3 REF 2003

### Legends

	Read instructions
	Expiration date
	LOT Charge number
	Keep in a dry place
	Store tightly closed keep
	Temperature range from 15 °C to 25 °C
	Manufacturer

CE 0029

### Warning

**H317**

May cause an allergic skin reaction.

**H319**

Causes serious eye irritation.

9999-2001-03

Revision issued

2015-07

**Dr. Ihde Dental AG**

Dorfplatz 11  
CH-8737 Gommiswald/SG  
Tel. +41 (0)55 293 23 23  
Fax +41 (0)55 293 23 00  
contact@implant.com  
www.implant.com

**Distribution by**

**Dr. Ihde Dental GmbH**

Erfurter Str. 19  
D-85386 Eching/München  
Tel. +49 (0)89 319 761-0  
Fax +49 (0)89 319 761-33  
info@ihde-dental.de  
www.ihde-dental.de

IHDENT® and GIZ® are registered trademarks.

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