Deutsch:

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English:

Instructions for use – Please read carefully

## IHDENT® GIZ® fil Type II

IHDENT® GIZ® fil Type II is a chemically curing, water-mixable glass-ionomer cement used as a filling material in restorative dentistry, as an orthodontic cement and as a cement for build-ups prosthodontics. IHDENT® GIZ® fil Type II is a flouride donor.

### Essential components

Powder: Aluminosilicate flouride, iron oxide pigment, polyacrylic acid

### Indications

- Class V restorations in cervical lesions or of mechanical / abrasive origin
- Class III restorations, especially in large lesions and in lower incisors
- Restorations of all classes in deciduous teeth
- Cementing of orthodontic appliances
- Cementing in prosthodontics
- Build-ups under crowns and bridges for crown preparations in prosthodontics
- Cavity linings of all kinds, especially for composite filling placed using the sandwich technique
- Carious root lesions
- As a fissure sealant for premolars and molars

### Application

#### Cavity preparation

Application in areas subject to high functional loads requires retentive preparations. Deep cavities above a vital dental pulp mandate the use of a calcium hydroxide liner. Class V cavity preparations must be at least 1 mm deep. Avoid knife-edge preparations.

1. Shake the powder vial well before dispensing.
2. Fill the water vial approximately halfway.
3. Remove powder from the powder vial using the measuring spoon and divide into three small heaps on the mixing block.
4. Drip water from the water vial: The prescribed mixing ratio of 1:1 (one drop: one measuring spoon) must be meticulously observed. Mix a minimum of 2 drops of water and 2 measuring spoons of powder. Using warm water accelerates the setting process.
5. Successively mix the three small heaps of powder into the liquid until the consistency of the filling material is homogeneous. Use an abrasion-resistant spatula (e.g. a plastic spatula) for this purpose.
6. Mix for a maximum of 30 seconds.
7. The filling material may be applied using any popular technique.
8. Once placed, fillings are sensitive to moisture. For this reason, apply a varnish (REF 2440) to the exposed filling surface.

Mixing time at 23 °C: 30 seconds  
 Working time: 2 minutes  
 Initial setting time: 4 - 5 minutes  
 Minimum time to finishing: 7 minutes from mixing

#### Finishing and polishing

IHDENT® GIZ® fil Type II fillings can be polished. Polishing can begin 7 minutes after mixing. Using small Arkansas stones or polishing discs at moderate speeds, using instruments of increasingly finer grit. Use petroleum jelly as a lubricant. Do not wet-polish until after 24 hours.

Once the polishing has been completed, the restoration can once again be protected with varnish (REF 2440) to protect the continuing hardening process.

#### Orthodontic use

Once the orthodontic band has been adapted, dry the appropriate tooth and carefully coat the band with IHDENT® GIZ® fil Type II on the inside. Then immediately place the band. You must remove any cement excess during the short viscous-plastic phase of the setting process. The position of the band must not be changed. To remove the band, loosen it and take the band off. Remove residual cement using a calculus removal unit or periodontal instruments.

#### Adverse effects

Glass-ionomer cements can be strongly irritating to the dental pulp. No such adverse effects have been reported for IHDENT® GIZ® fil Type II to date.

#### Contraindications

The use of IHDENT® GIZ® fil Type II is contraindicated for direct application to dentin close to the dental pulp or to the exposed dental pulp itself, for restorations subject to high masticatory loads and in patients known to be allergic to one to the ingredients of IHDENT® GIZ® fil Type II.

#### Interactions

None known.

#### Warnings

When glass-ionomer powder or mixed cement come in contact with the eye, chemical burns may be caused by the polyacrylic acid. If contact with the eye has occurred, immediately and carefully rinse with copious amounts of water and consult a knowledgeable ophthalmologist.

#### Notes

Rinse plastic or metal spatulas with water after use and remove any material residue. The powder vial contains a globe that extracts residual liquid from the powder. Do not remove or open the globe.

**Keep out of the reach of children. Intended for use by dentists only.**

### Package sizes + Accessories

Package, 50 g powder	
Shade A2 universal	REF 2051
A3 universal	REF 2053
A3.5 brown	REF 2057

### Starter pack

5 shades, 10 g each, varnish, accessories REF 2010

### Sample pack, 10 g powder

Shade A2 universal REF 2011

### Accessories

Mixing block	REF 2070
Measuring spoon	REF 2020
Water vial	REF 2025
Varnish 30 ml	REF 2440
Conditioner 30 ml	REF 2330

### Storage / Shelf life

Store tightly closed in a cool and dry place. Maximum storage temperature: 25 °C. Check the bar code on the packaging or on the vial label for batch number and expiry date.

### Legends

-  Consult instructions for use
-  Expiration date
-  (Gamma) -sterilized
-  Only use once
-  LOT Charge number
-  Store protected against the sun (Conditioner)
-  Keep in a dry place
-  Store tightly closed keep
-  Temperatures up to 25 °C
-  Manufacturer

# CE 0029

ISO 9917:2007

9999-2010-03

Revision issued 2015-08

 **Dr. Ihde Dental AG**  
 Dorfplatz 11  
 CH-8737 Gommiswald/SG  
 Tel. +41 (0)55 293 23 23  
 Fax +41 (0)55 293 23 00  
 contact@implant.com  
 www.implant.com

**Distribution by Dr. Ihde Dental GmbH**  
 Erfurter Str. 19  
 D-85386 Eching/München  
 Tel. +49 (0)89 319 761-0  
 Fax +49 (0)89 319 761-33  
 info@ihde-dental.de  
 www.ihde-dental.de

IHDENT® and GIZ® are registered trademarks.

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