Deutsch:

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English:

Instructions for use – Please read carefully

IHDENT® GIZ® fil XN

IHDENT® GIZ® fil XN is a light cured reinforced glass-ionomer filling cement. Beside its high fluoride release and the excellent biocompatibility as glass-ionomer cement IHDENT® GIZ® fil XN has also good chemical bonding to dentine and enamel and a tight seal at the dentinal margins. Because of its radiopacity it ensures easy postoperative diagnosis.

Essential components

Mixture of strontium-fluoro-alumino silicate, glasspowder and pigments.

Application areas

- uniplanar fillings of class I, III and V
- primary tooth fillings
- linings for composite fillings
- core build-ups
- small class II fillings

Application

Preparing the tooth

Grind the tooth using standard techniques. Do not prepare thin edges. If desired, place a matrix band. Areas close to the pulp should be covered with a small amount of a hard setting calcium hydroxide liner (REF 8768).

As an option (in order to increase the bond between filling and tooth surface) a conditioner may be applied for 20 seconds (Conditioner REF 2330). For optimal chemical bonding IHDENT® GIZ® fil XN can be applied for 20 seconds to the bonding surfaces to remove the preparation smear layer. Rinse IHDENT® GIZ® fil XN thoroughly with water and dry gently. Do not desiccate.

Dispensing and Mixing

The powder / liquid ratio to achieve a suitable consistency is 2.6 / 1.0. This can be obtained by mixing 1 level scoop of powder and 1 drop of liquid.

For accurate dispensing of IHDENT® GIZ® fil XN Powder shake the bottle to loosen the powder. Overfill the spoon with the powder, level the powder by using the collar at the top of the bottle and carry it onto the mixing pad provided. Avoid compressing powder into the spoon with the inside wall of the bottle. For dispensing of IHDENT® GIZ® fil XN Liquid turn the bottle vertically with the tip about 5 cm above the mixing pad. Steady your hand and squeeze the bottle gently to dispense the drops. If any bubbles are present, lightly tap the bottle with the fingers holding it. Discard drops that contain bubbles and are obviously not full-sized. Use a small spatula to rapidly mix half the powder into the liquid. Then mix the remaining powder into the mixed paste to achieve a glossy appearance. Total mixing time is 30 seconds. After use, close both powder and liquid bottles.

Filling

Mix the required amount of IHDENT® GIZ® fil XN and apply it into the prepared cavity by using a suitable placement instrument. Please see to it, that no air bubbles will be incorporated. IHDENT® GIZ® fil XN should be placed in the cavity within the working time (ca. 1:30 minutes from start of mixing at 23 °C). If desired, a matrix band may be used to form the contour. Net setting time is about 4 minutes.

Note: Higher temperatures will shorten the working time; lower temperatures will prolong the working time.

An overextended working time will cause the loss of adhesion to the dental enamel and the dentine. Remove the matrix when the cement has achieved clinical set (approx. 4 minutes after application).

Finishing

Final finishing and polishing can begin from about 6 minutes after start of mixing. Smooth with fine diamond burs, then polish with finishing and polishing discs with graded grain size.

Optional apply a thin layer of Zero Glas Coat to the final finished surface of the restoration. Immediately light cure for 20 seconds.

Close bottle immediately after use. Instruct the patient not to expose the restoration to any pressure for one hour.

Conclusive Notes

The products are to be applied only by a dental professional in the manner as described in this instruction. Do not use the products with patients who show an allergy to the material. In case of allergic reactions immediately stop the application, and advise the patient to consult a physician.

Do not allow the liquid or the mixture to contact the oral tissues or skin. In case of contact, remove the material with absorbent cotton soaked in alcohol and rinse with water. Avoid eye contact of the mixture. In case of contact, immediately flush with water and seek medical treatment.

Do not mix the powder or liquid of IHDENT® GIZ® fil XN with any other restorative cement components.

Precautionary statements

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P262 Do not get in eyes, on skin, or on clothing.

Storage / Shelf life

Store the products at room temperature. Temperature should not exceed 25 °C. Do not use after expire date. Refer to the barcode on the packaging and syringe label for Batch Number and Expiry Date.

Package sizes


15 g powder and 8 ml liquid

Accessories: 1 dosage spoon, 1 mixing block, 1 mixing plate

Shade: Universal **A2**

REF 2082

Legend

	Read instructions
	Expiry date
	LOT Charge number
	Only use once
	Keep in a dry place
	Store tightly closed keep
	Temperature range from 4 °C to 25 °C
	Manufacturer

CE 0029

Warning

H315

Causes skin irritation.

H319

Causes serious eye irritation.

H335

May cause respiratory irritation.

9999-2082-04

Revision issued 2015-08

Dr. Ihde Dental AG

Dorfplatz 11
CH-8737 Gommiswald/SG
Tel. +41 (0)55 293 23 23
Fax +41 (0)55 293 23 00
contact@implant.com
www.implant.com

Distribution by

Dr. Ihde Dental GmbH

Erfurter Str. 19
D-85386 Eching/München
Tel. +49 (0)89 319 761-0
Fax +49 (0)89 319 761-33
info@ihde-dental.de
www.ihde-dental.de

IHDENT® and GIZ® are registered trademarks.

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