Deutsch:

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English:

Instructions for use – Please read carefully

IHDENT® TempBase® Perfekt

IHDENT® TempBase® Perfekt is an automatically mixed two component material in cartridges based on multifunctional methacrylic esters.

Essential components

Urethane dimethacrylate, aromatic dimethacrylate, Glycolmethacrylat.

Application areas

IHDENT® TempBase® Perfekt is a very easy and comfortable material for creating temporary crowns or bridges, onlays and veneers. Because of its flexibility, the material is especially suitable for longer bridge spans. IHDENT® TempBase® Perfekt shows fluorescence in UV-light.

IHDENT® TempBase® Perfekt is free of methyl methacrylate. Its temperature derived while curing is lower than 40 °C. As a temporary crown or bridge it protects the prepared teeth against external influences and preserves the occlusion. It shows increased transverse strength, good abrasion resistance and low polymerization shrinkage. It fits perfect. Showing good polishability, good color stability and fluorescence it looks perfect. Temporaries made from IHDENT® TempBase® Perfekt is easy to repair. Cured material can be repaired with any chemicaly or light curing composite building up again a chemical link. In nearly every property IHDENT® TempBase® Perfekt is superior to the common materials based on monomer/polymer systems. It allows thin edges and production of longer bridge spans.

Impression taking

Before preparing of stumps or before a planned extraction of a tooth a situation impression is made by addition curing silicones (long storage stability!) or by alginates. Interdental areas are carved out, in molar areas with teeth absent it may be necessary to cut a groove in the impression between the abutments to create a bridge-like connection between the tooth units.

Preparing of the cartridge

The turn cap of the cartridge is taken off and thrown away (do not use it again!). It is substituted by the supplied mixing cannula by attaching it with a sideways turning movement, locking the cannula in place.

The application gun is loaded with the prepared cartridge and is ready for application. The initial extrusion from the mixing cannula (about the size of a pea) should be discarded. Then the following mix will be perfect. This must be done for each new mix. Leave the used mixing cannula on the cartridge. It serves as a cap.

Application

IHDENT® TempBase® Perfekt is automatically mixed when dispensed with slight and even pressure directly into the situation impression made before. Filling should occur from bottom upward to prevent voids.

Note

- Load the situation impression with IHDENT® TempBase® Perfekt
- Seat the impression onto the prepared areas of the teeth. Apply in a maximum of 60 seconds (working time at 23 °C).
- After 1 - 2 minutes (setting time in mouth) the material shows a semi-hardened but still elastic condition and can be removed from the teeth (alternatively: 3 - 4 minutes setting time on the model)

Note:

The temporary crown and bridge substitute can only be removed without destruction during the elastic state.

Processing

After removal of the temporary crown or bridge from the impression (or from the tooth preparation) excess material and proximal undercuts are removed. Preferably, the temporary is then post cured in warm (45 °C - 55 °C) water (e.g. in a hot cure polymerization device). If this is not possible or not desired, after 6 minutes from beginning of application the temporary acrylic can be worked out with rotary instruments and can be polished to high gloss.

Do not breathe polishing dust; use suitable mouth protective device or aspiration!

Remark

The oxygen inhibited smear layer on the surface caused by the acrylic system can not be fully avoided and should be removed before working out or cementing. It can easily be removed by alcohol or other suitable solvents.

Cementing / Repairs of the temporary

IHDENT® TempBase® Perfekt-Provisorien temporary crowns or bridges preferably should be cemented with an eugenol-free temporary cement (e. g. JetBond II REF 2126). In case of using eugenol-containing cements it should be considered that later used acrylate based cements could be hindered in hardening (inhibition of curing reaction by traces of remained eugenol).

Repairs

IHDENT® TempBase® Perfekt temporaries show high mechanical strength. However, if a temporary breaks the following procedures are recommended:

Fracture shortly after production:

Both ends of the fracture are refixed with freshly extruded IHDENT® TempBase® Perfekt. Enamelite® REF 4050 can be applied alternatively.

Fracture of longer existing temporary:

The areas of fracture are cleaned and roughened and provided with some mechanical retentions. Thus prepared, the fractured areas can be joint with freshly extruded IHDENT® TempBase® Perfekt. To enhance the curing place for some minutes in warm water. A flowable light cure composite (e.g. Enamelite® REF 4050) can be applied alternatively. To ensure a complete cure, do not apply Enamelite® in layers exceeding 2 mm. Light cure each layer for 40 seconds.

Important note

Non hardened IHDENT® TempBase® Perfekt can easily be removed with alcohol or other suitable solvents.

Warnings

IHDENT® TempBase® Perfekt is free of methyl methacrylate but contains other methacrylates.

With susceptible patients, sensitization to IHDENT® TempBase® Perfekt cannot be excluded. IHDENT® TempBase® Perfekt should not be used any more, if allergic reactions are observed. Do not use for patients with allergic reactions against acrylates.

Avoid contact with skin, mucous membrane and eyes. If the material comes into contact with skin, immediately wash with water and soap. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.

Keep out of the reach of children.

Intended for use by dentists only.

Storage

All times mentioned here refer to application at 25 °C. Check the bar code for batch number and expiry date. Do not use after expiry date.

Presentation

Regular package

1 cartridge 50 ml / shade universal A2 / 15 cannulas
mixing ratio **4:1** REF 7045

Refill pack

3 cartridges à 50 ml / shade universal A2
mixing ratio **4:1** REF 7046

Accessories

mixing cannulas blue pack of 25 REF 7062
safety mixing gun **4:1** REF 7747

Legends

	Read instructions
	Expiration date
	LOT Charge number
	Keep in a dry place
	Store tightly closed keep
	Temperatures up to 25 °C
	Manufacturer

CE 0029

9999-7045-03

Revision issued 2015-07

Dr. Ihde Dental AG

Dorfplatz 11
CH-8737 Gommiswald/SG
Tel. +41 (0)55 293 23 23
Fax +41 (0)55 293 23 00
contact@implant.com
www.implant.com

Distribution by

Dr. Ihde Dental GmbH

Erfurter Str. 19
D-85386 Eching/München
Tel. +49 (0)89 319 761-0
Fax +49 (0)89 319 761-33
info@ihde-dental.de
www.ihde-dental.de

IHDENT®, TempBase® and Enamelite® are registered trademarks.

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