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Manufacturer:

P.L. Superior Dental Materials GmbH, 22607 Hamburg, Germany

Dual cure flowable core build-up and post cementing composite with nano-technology

HARDCORE® DC is a dual cure high radiopaque flowable microhybrid composite material with nano particles (nano-zirconium dioxide, nano-calciumfluoride) for core-build-up and cementation of posts. Due to its excellent mechanical properties final crown preparation can be carried out more precise. **HARDCORE® DC** is also suitable for cementing of crowns and bridges, inlays and onlays if there are no high esthetic requirements.

Delivered in 1:1 AUTOMIX-cartridges or MINIMIX-syringes **HARDCORE® DC** can be easily dispensed and applied directly. **HARDCORE® DC** exhibits a short setting time without high heat generation. The dual cure properties enables the dentist, also to carry out cementations and core build-ups in cases where a light cure cannot be guaranteed to be sufficient.

HARDCORE® DC meets the requirements of **DIN EN ISO 4049**.

Side Effects

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection the cavity floor in cases of deep excavations should be covered with a thin layer of calciumhydroxide material.

Incompatibility with Other Materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite material in proximity of eugenol containing products, nor let the composite allow coming into contact with materials containing eugenol.

Preparing of the cartridge (HARDCORE® DC)

Take off the turn cap of the AUTOMIX-cartridge and throw it away (**do not use it again!**). It is substituted by the special 1:1 mixing cannula by attaching it with a sideways turning movement, locking the cannula in place.

The application gun is easy to load with the prepared cartridge and is ready for application.

Note:

The initial extrusion from the mixing cannula (about the size of a pea) should be discarded. Then the following mix will be perfect. This must be done for each new mix.

Preparing the MINIMIX-Syringe (HARDCORE® DC XPress)

Remove the cap of the MINIMIX-syringe and throw it away (**do not use it again!**). It is replaced by a special 1:1 mixing cannula. Turn the cannula 90° until it locks in position. The Material is now ready for application.

Note:

Discard the first 2-3 mm of the extruded material (about the size of a peppercorn). This must be done for each new mix.

Store used AUTOMIX-cartridge/MINIMIX-syringe with fixed used mixing cannula in the dark.

The working time (23°C (74°F)) of **HARDCORE® DC** in the self cure mode is **1:30 minutes** from start of mixing.

1. Post Cementation**1.1. Isolation**

Use of a rubber dam to isolate the tooth is strongly recommended.

1.2. Root Canal Preparation

Refer to directions of the selected post manufacturer.

Remove all existing old restorations and decay from the tooth. Prepare and clean the root canal with e.g. sodium hypochlorite solution, rinse and remove excess solution from the canal with a soft paper tip.

Apply a suitable dual cure adhesive for enamel and dentine bonding according to the corresponding manufacturer instructions.

The application of the adhesive Dentamed® E&B® 2 is recommended. Follow the instructions for Dentamed® E&B® 2.

1.3. Post Cementation

Prepare the selected post according to manufacturer directions.

HARDCORE® DC is applied into the prepared root canal and onto the post. Seat the post careful into the canal and maintain firm pressure until the post is seated. **HARDCORE® DC** self-cures within **3:30 minutes**. For post stabilization light cure the coronal part of the cemented post for **20 seconds**.

As soon as the **HARDCORE® DC** has set proceed with the core-build-up procedure.

2. Core-Build-Up**2.1 Isolation**

Use of a rubber dam to isolate the tooth is strongly recommended.

2.2. Cavity Preparation

Remove all existing old restorations and decay from the tooth. If necessary place any pins or posts. Refer to directions of the selected post manufacturer.

2.3. Pulp Protection

Cavity floor of deep excavations should be covered with a thin layer of calciumhydroxide material.

2.4. Application of a Bonding Agent

Apply a suitable dual cure adhesive for enamel and dentine bonding according to the corresponding manufacturer instructions.

The application of the adhesive Dentamed® E&B® 2 is recommended. Follow the instructions for Dentamed® E&B® 2.

Note:

It is essential that the primed dentine and enamel surfaces are dry and contaminant free for the application of **HARDCORE® DC**.

2.5. Application

Place the mixing cannula directly into the preparation and press out the paste.

HARDCORE® DC is automatically mixed when dispensed with slight and even pressure. Filling should occur from bottom upwards to prevent voids. To facilitate placement of **HARDCORE® DC** place a matrix band around the prepared tooth.

HARDCORE® DC may be contoured by using a composite instrument. A flat-ended interproximal carver is recommended.

Place **HARDCORE® DC** directly into the preparation and allow the system to self cure for **3:30 minutes**. After that the material should be light-cured for **40 seconds**. With this technique an optimum of physical properties will be obtained.

An explorer can be used to test that the **HARDCORE® DC** has completely set. Remove the matrix not earlier than the material has set.

Final crown preparation on the **HARDCORE® DC** can be carried out by using crown preparation burs.

Additional Notes

- Do not use any resin to adjust viscosity of composite restorative material.
- Do not use in combination with substances containing eugenol because eugenol impairs the hardening of the composite. Do neither store the composite material in proximity of eugenol containing products, nor let the composite allow to come into contact with materials containing eugenol.
- Contact of resin pastes with skin and gingival tissue should be avoided, especially by anyone having known resin allergies.

Storage

Do not store above 20 °C (68 °F). Store unopened material in the refrigerator.

Opened cartridges have to be used up within 3 months.

Do not use after expiry date.

Warranty

P.L. Superior Dental Materials GmbH warrants this product will be free from defects in material and manufacture. P.L. Superior Dental Materials GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and P.L. Superior Dental Materials GmbH's sole obligation shall be repair or replacement of the P.L. Superior Dental Materials GmbH product.

Limitation of Liability

Except where prohibited by law, P.L. Superior Dental Materials GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep away from children!

For dental use only!

Distributor:

Dr. Ihde Dental GmbH
Erfurter Straße 19
D-85386 Eching
Tel. +49 89 319 761-0
info@ihde-dental.de